




TRansition paths to sUustainable  
legume-based systems in EUrope

## Antioxidant assay 1: Ascorbate peroxidase (APX)

---

**Date:** September, 30<sup>th</sup> 2018



This Project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 727973

[www.true-project.eu](http://www.true-project.eu)



---

## Content

<b>1 Summary Information .....</b>	<b>3</b>
1.1 Partner Summary.....	3
1.2 SOP Summary .....	4
<b>2 Protocol Steps .....</b>	<b>5</b>
<b>3 Linked SOPs .....</b>	<b>10</b>
<b>4 Disclaimer .....</b>	<b>10</b>
<b>5 Copyright.....</b>	<b>10</b>
<b>6 Citation .....</b>	<b>10</b>



## 1 Summary Information

### 1.1 Partner Summary

SOP Code	EU_TRUE_SOP_010
TRUE Partner Acronym	AUA
Primary Author	Ntatsi, Georgia ( <a href="mailto:ntatsi@aua.gr">ntatsi@aua.gr</a> )
Other Authors	Savvas, Dimitrios
Linked Reference and Hyperlink (if available)	Nakano, Y., Asada, K., 1981. Hydrogen peroxide is scavenged by ascorbate-specific peroxidase in spinach chloroplasts. Plant and Cell Physiology, Volume 22, Issue 5, August 1981, Pages 867-880 <a href="https://doi.org/10.1093/oxfordjournals.pcp.a076232">https://doi.org/10.1093/oxfordjournals.pcp.a076232</a>
Associated files to use with the SOP [and function]	NA



---

## 1.2 SOP Summary

### Title

**Antioxidant assay 1: Ascorbate peroxidase (APX)**

### Brief description

Ascorbate (AsA) is a major antioxidant and free-radical scavenger in plants and is considered to be of paramount importance as an electron donor for  $H_2O_2$  detoxifications *via* ascorbate peroxidase (APX) in plant cells (Ahmad *et al.*, 2010)

## 2 Protocol Steps

### 1<sup>st</sup> stage: Enzyme Extraction

0.1 g of each pulverized, frozen sample (leaf and root samples) was homogenised with ice-cold **25 mM HEPES buffer (pH 7.8)**. The HEPES buffer contained **0.2 mM EDTA, 2 mM ascorbate and 2 % (w/v) polyvinylpyrrolidon (PVP)**. The homogenates were centrifuged at 4 °C and 14,000 rpm. The supernatants obtained were used for enzyme analysis. The procedure for enzyme extraction was carried out at 0-4 °C.

### Reagents

To prepare **1 L** of the necessary reagents, the following equation was used:

$$m = C * V * M \quad (\text{equation 1})$$

with m = weight (g)  
C = molarity (mol/L)  
V = volume to be prepared (L)  
MB = molecular weight (g/mol)

- **25 mM HEPES buffer (pH 7.8)**

Weigh 5.9575 g of HEPES (MB = 238.30) in a 1 L volumetric flask and adjust to 1 L with H<sub>2</sub>O ( $m = 0.025 * 1 * 238.30$ ). Adjust the pH to 7.8 using either an HCl or NaOH solution.

- **0.2 mM EDTA**

Weigh 0.0745 g of EDTA (MB = 372.24) in a 1 L volumetric flask and adjust to 1 L with H<sub>2</sub>O ( $m = 0.0002 * 1 * 372.24$ ).

- **2 mM ascorbate (AsA)**

Weigh 0.352 g of Ascorbate (MB = 176.12) in a 1 L volumetric flask and adjust to 1 L with H<sub>2</sub>O ( $m = 0.002 * 1 * 176.12$ ).

- **2% (W/V) polyvinylpyrrolidon PVP**

- 2 g in 100 mL

## 2<sup>nd</sup> stage: Measuring Ascorbate peroxidase

The ascorbate peroxidase (APX) activity was measured in the leaf-samples in accordance to the method of Nakano and Asada (1981), with a few modifications due to ascorbate oxidation and therefore the decrease in absorbance at 290 nm .

200  $\mu$ L of reaction mixture containing **25 mM potassium phosphate buffer (pH 7.0), 0.5 mM ascorbate, 0.1 mM EDTA and 0.1 mM hydrogen peroxide** and 10  $\mu$ L plant extract were used. The reaction was started by adding the hydrogen peroxide and the absorbance at 290 nm was measured over five minutes (7-14 minutes better) with the Synergy HT 96- position microplate spectrophotometer (Biotek, Winooski, USA). The absorbance of the mixture before hydrogen peroxide was added served as blank. To transfer the data in an excel format, the programme KC4 was used. For the calculation of APX-activity the absorbance coefficient of reduced ascorbate ( $E = 2.8 \text{ mM} \cdot \text{cm}^{-1}$ ) was used (Nakano and Asada, 1981).

### Reagents for the preparation of the reaction mixture:

- **25 mM potassium phosphate buffer (pH 7.0)**

#### **1M $\text{K}_2\text{HPO}_4$**

Since the molar mass of  $\text{K}_2\text{HPO}_4$  is 174.1760 g/mol, a 1 M solution would be 174.1760 g in 1L.

#### **1M $\text{KH}_2\text{PO}_4$**

Since the molar mass of  $\text{KH}_2\text{PO}_4$  is 136.0857 g/mol, a 1 M solution would be 136.0857 g in 1L.

### Preparation of 1L of 0.1 M = 100 mM potassium phosphate buffer at 25 °C

Desired pH	Volume of 1M $\text{K}_2\text{HPO}_4$ (mL)	Volume of 1M $\text{KH}_2\text{PO}_4$ (mL)
7.0	61.5	38.5

Hence, for the preparation of 100 mL of 0.025 M = 25 mM potassium phosphate buffer at 25 °C

Desired pH	Volume of 1M $\text{K}_2\text{HPO}_4$ (mL)	Volume of 1M $\text{KH}_2\text{PO}_4$ (mL)
7.0	1.5375	0.9625

- **0.5 mM ascorbate**

To prepare 100 mL of a 0.5 mM solution of ascorbate, use 8.806 mg of Ascorbate ( $M_B = 176.12$ ) ( $m = 0.0005 \cdot 100 \cdot 176.12$  from equation 1).

- **0.1 mM  $\text{H}_2\text{O}_2$**

First determine the molarity (C) of your commercial  $\text{H}_2\text{O}_2$  solution (mol/g).

Using the equation 2 below:

$$C = (m/MB) * V \quad (\text{equation 2})$$

with  $m$  = weight (g) =  $D * \% \text{ solution}$   
 $MB$  = Molecular weight (g/mol)  
 $V$  = Volume (L)

The molarity of a 30% solution of  $H_2O_2$  can be calculated using:

$$C = (D * \% \text{ solution}) / MB * 1,000 \quad (\text{equation 3})$$

$MB = 34.01$   
 $\% \text{ of solution} = 30$   
 $\text{Density} = 1.1 \text{ g/cm}^3 = 1,100 \text{ g/L}$

Hence the molarity of 30 % solution of  $H_2O_2$  is 9.79 M or 9,790 mM ( $1.1 * 30/100$ ) /  $34.01 * 1,000$ ).

To prepare 1 L of a 0.1 mM solution, use the following equation:

$$C_1V_1 = C_2V_2$$

with  $C_1 = 9,790 \text{ mM}$   
 $C_2 = 0.1 \text{ mM}$   
 $V_2 = 1,000 \text{ mL}$

$$V_1 = 0.1 * 1,000 / 9,790 = 0.0102 \text{ mL}$$

Similarly, to prepare a 100 mL of a 10 mM solution, use 0.102 mL of your  $H_2O_2$  commercial solution ( $V_1 = 10 * 100 / 9,790$ )

Then prepare a diluted solution of 0.1 mM (1 mL of 10 mM into 100 mL  $H_2O$ )  
<http://www.graphpad.com/quickcalcs/molarityform.cfm>

- **0.1 mM EDTA**

To prepare 100 mL of a 0.1 mM solution of EDTA, use 3.7224 mg of EDTA ( $MB = 372.24$ ) ( $m = 0.0001 * 100 * 372.24$  from equation 1 above).

### Preparation of 100 mL of the reaction mixture

Solution	Stock Solution Molarity	Quantity of stock
25 mM $K_2HPO_4$	1 M $K_2HPO_4$	1.5375 mL
25 mM $KH_2PO_4$	1 M $KH_2PO_4$	0.9625 mL
0.5 mM ascorbate		8.806 mg
0.1 mM $H_2O_2$	10 mM	1 mL

0.1 mM EDTA		3.7224 mg
H <sub>2</sub> O		96.5 mL
<b>Total Volume</b>		<b>100 mL</b>

### 3<sup>rd</sup> stage: Calculations

First with the absorbance measured, calculate the C using equation 4:

$$C = \Delta A / (\epsilon \cdot d) = \Delta A / [2.8 \cdot 10^3 \cdot d] \quad (\text{equation 4})$$

with C = molarity of ascorbate (mol/L)

$\Delta A$  = Absorption Difference

e = Extinction absorbance coefficient of reduced ascorbate:

$$\epsilon = 2.8 \text{ mM} \cdot \text{cm}^{-1} = 2.8 \cdot 10^3 \text{ L/mol} \cdot \text{cm}$$

D = Light path (cm) (0.1 cm plastic plate, 0.2 cm glass plate)

To calculate the ascorbate peroxidase activity, use equation 5 below:

$$U \text{ (mol/g/min)} = (V \cdot C) / (\Delta t \cdot FM) \quad (\text{equation 5})$$

with U = Activity (mol min<sup>-1</sup> gr<sup>-1</sup>)

V = Assay volume (L)

$\Delta t = t_1 - t_0$  = (time between measurements) (min)

FM = weight of the fresh material (g)

Or by substituting C in equation 5 for equation 4:

$$U = (V \cdot C) / (\Delta t \cdot FM) = [V \text{ (L)} \cdot C] / [\Delta t \cdot FM]$$

From this, you can also calculate the ascorbate peroxidase activity *per* protein (U') by dividing U by the protein content in grams (Equation 6).

$$U' \text{ (mol gr}^{-1} \text{ protein min}^{-1}) = U / \text{Protein content (g)}^1 \quad (\text{equation 6})$$

<sup>1</sup> See EU\_TRUE\_SOP\_015 on how to determine Protein Content.





---

## Reference

Parvaiz Ahmad, Cheruth Abdul Jaleel, Mohamed A. Salem, Gowher Nabi & Satyawati Sharma (2010)  
Roles of enzymatic and nonenzymatic antioxidants in plants during abiotic stress, Critical Reviews  
in Biotechnology, 30:3, 161-175, DOI: 10.3109/07388550903524243

### 3 Linked SOPs

SOP Code	SOP Function
EU_TRUE_SOP_011	Antioxidant assay 2: Catalase (CAT)
EU_TRUE_SOP_012	Antioxidant assay 3: Guaiacol peroxidase (G-POD)
EU_TRUE_SOP_013	Antioxidant assay 4: Glutathione reductase (GR)
EU_TRUE_SOP_014	Antioxidant assay 5: Superoxide dismutase (SOD)
EU_TRUE_SOP_015	Bradford protein assay

### 4 Disclaimer

The information presented here has been thoroughly researched and is believed to be accurate and correct. However, the authors cannot be held legally responsible for any errors. There are no warranties, expressed or implied, made with respect to the information provided. The authors will not be liable for any direct, indirect, special, incidental or consequential damages arising out of the use or inability to use the content of this publication.

### 5 Copyright

© All rights reserved. Reproduction and dissemination of material presented here for research, educational or other non-commercial purposes are authorised without any prior written permission from the copyright holders provided the source is fully acknowledged. Reproduction of material for sale or other commercial purposes is prohibited.

### 6 Citation

Please cite this report as follows:

Ntatsi, G., Savvas, D. (2018) Standard Operating Procedure 06: Antioxidant assay 4: Ascorbate peroxidase (APX). Developed by the EU-H2020 project TRUE ('Transition paths to sustainable legume-based systems in Europe'), funded by the European Union's Horizon 2020 Research and Innovation programme under Grant Agreement Number 727973. Available online at: [www.true-project.eu](http://www.true-project.eu).