




TRansition paths to sUustainable  
legume-based systems in EEurope

## Soil Organic Carbon (OC)

---

**Date:** September, 30<sup>th</sup> 2018



This Project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 727973

[www.true-project.eu](http://www.true-project.eu)



---

## Content

<b>1 Summary Information</b>	<b>3</b>
1.1 Partner Summary	3
1.2 SOP Summary	4
<b>2 Protocol Steps</b>	<b>5</b>
Spectrophotometric determination of Phosphorus content in plant tissues (MURPHY-RILEY)	错误!未定义书签。
Determination of K in plant tissues	错误!未定义书签。
<b>3 Linked SOPs</b>	<b>9</b>
<b>4 Disclaimer</b>	<b>9</b>
<b>5 Copyright</b>	<b>9</b>
<b>6 Citation</b>	<b>9</b>



## 1 Summary Information

### 1.1 Partner Summary

<b>SOP Code</b>	EU_TRUE_SOP_026
<b>TRUE Partner Acronym</b>	AUA
<b>Primary Author</b>	Savvas, Dimitrios ( <a href="mailto:dsavvas@aua.gr">dsavvas@aua.gr</a> )
<b>Other Authors</b>	Karavidas I., Ntatsi G.
<b>Linked Reference and Hyperlink (if available)</b>	Walkley, A. and I.A. Black. 1934. An examination of the Degtjareff method for determining soil organic matter, and a proposed modification of the chromic acid titration method. Soil Sci.37:29-38 <a href="http://dx.doi.org/10.1097/00010694-193401000-00003">http://dx.doi.org/10.1097/00010694-193401000-00003</a>
<b>Associated files to use with the SOP [and function]</b>	<i>Not applicable</i>

## 1.2 SOP Summary

### Title

#### Soil Organic Carbon (OC)

### Brief description

Soil organic carbon improve soil structure, water retention, soil biodiversity, buffering capacity and the cycling and storage of plant nutrients. In addition, the mineralization of soil organic carbon provides plant nutrients, especially nitrogen, phosphorus and sulfur (S), as well as with micronutrients. Therefore, the measure of soil organic carbon levels is necessary to determine the soil fertility. The Walkley-Black method was chosen as it is a simplified method, which provides accurate measurements on soils with high organic carbon content.

#### Application and Principle

The Walkley-Black (WB) method (Walkley and Black, 1934) determines soil organic carbon (OC) by oxidising carbon with acidic dichromate ( $K_2Cr_2O_7$ ). The oxidation step is followed by titration of excess dichromate with ferrous sulfate. The OC is calculated from the difference between the total dichromate added and the amount of dichromate left unreacted after the OC oxidation. The method provides an estimate of soil Organic Matter from OC determination assuming that the fraction of soil OC reacted and the amount of OC in soil OM. This method was widely considered a standard method for OM after its development, but the advance of thermal decomposition instruments has resulted in its limited use. Instead loss on ignition and thermal decomposition methods are more widely used due to their ability to analyse samples in a more time efficient manner.

---

## 2 Protocol Steps

### Equipment and Apparatus

1. Analytical balance with 0.01 g resolution
2. 250-mL wide mouth graduated Erlenmeyer flasks
3. Repipetter to deliver 10 mL volume
4. Fume Hood
5. Titration stand and burette
6. Stir plate with light
7. Stirring rods
8. Weighing vessel

### Reagents

1. **Potassium dichromate ( $K_2Cr_2O_7$ )**, 0.167 M, 1 N  
Dissolve 98.08 g of oven-dried potassium dichromate in approximately 1500 mL of deionized water and dilute to 2 L. After preparation, transfer to a clean glass bottle for use with a repipetter to deliver 10 mL. Do not mix old potassium dichromate solution with the new solution.
2. **Ferrous sulfate ( $FeSO_4 \cdot 7H_2O$ )**, 0.5 M, 0.5 N  
Dissolve 278.02 g of ferrous sulfate in approximately 1500 mL of deionized water. Carefully add 30 mL of concentrated sulfuric acid, mix, cool, and dilute to 2 L. After preparation, this solution may be transferred to a clean 8-L plastic carboy. Do not mix old ferrous sulfate solution with the new solution. The tubing, stopcock, and attachments to the burette should be rinsed three times with new ferrous sulfate solution before titrating any blanks or samples. Prepare a new solution every 30 days.
3. **Ferroine** (1,10-Phenanthroline-ferrous sulfate-complex) indicator solution, 0.025 M  
Available from chemical supply companies.
4. **Sulfuric acid ( $H_2SO_4$ )**  
Concentrated and not less than 96%.

### Procedure

1. Titrate blank samples without soil before proceeding with any unknown samples. Add 10 mL of potassium dichromate into a 250 mL Erlenmeyer flask. Titrate with ferrous sulfate. Titrant volume should be approximately 20 mL. Record the exact volume of titrant to the nearest 0.1 mL. If the titrant volume of the two blanks is not within 0.2 mL, clean the burette and associated tubing. Reanalyse two more blanks.
2. Weigh 1.00g of processed mineral soil (dried, <2 mm) into a 250-mL wide mouth graduated Erlenmeyer flask. If soil weight is not exactly 1.00 g, record weights for later calculation.

3. Pipette 10 mL of potassium dichromate into each flask and mix carefully by rotating the flask to wet all of the soil.
4. Under a fume hood, carefully add 20 mL of sulfuric acid to each flask and mix gently.
5. Allow flasks to stand for 5 min under the fume hood.
6. Add deionized water to each flask so the final volume is approximately 125 mL. Mix by swirling gently.
7. Allow the samples to cool and return to room temperature. If volume is less than 125 mL after 30 min, add more deionized water to achieve 125 mL volume.
8. Add 5 or 6 drops of ferroine indicator solution and immediately titrate with ferrous sulfate. Use a mixing bar to thoroughly stir the sample as it is titrated. As the titration proceeds, the solution will take on a green colour that will change abruptly to reddish-brown when the endpoint of the titration is reached.
9. Record each volumetric reading to the nearest 0.1 mL.

### Calculations

1. Soil organic C (OC) is determined from the  $m_{eq}$  of  $K_2Cr_2O_7$  (in mL) reacting with the soil to oxidize organic C, which is the difference between total  $m_{eq}$  of  $K_2Cr_2O_7$  (in mL) added to soil and the  $m_{eq}$  of  $FeSO_4$  that titrated the remaining  $K_2Cr_2O_7$  after the reaction.

$$OC (\%) = (m_{eq} K_2Cr_2O_7 - m_{eq} FeSO_4) \times [0.003 (g.C.m_{eq}^{-1})] \times cf \div g \text{ soil} \times 100$$

The oxidation state of C in organic matter is assumed to be 0. The C in organic matter is oxidised to  $CO_2$  with an oxidation state of +4. Thus, the equivalent weight of C is 3 g ( $12 \text{ g C} \div \text{mole} \times \text{mole C} \div 4 \text{ eq C} = 3 \text{ g C}_{eq}^{-1} = 0.003 \text{ g C } m_{eq}^{-1}$ ). An oxidation correction factor (cf) is required since not all the organic C in soil is oxidized with room temperature oxidation. Correction factors can range from 1.14 to 1.32.

2. The  $m_{eq}$  of  $K_2Cr_2O_7$  reacting with soil to oxidise OC is determined from the following equation.

$$(m_{eq} K_2Cr_2O_7 - m_{eq} FeSO_4) = (mLB - mLS) \times 0.5 N$$

with **mLB** = average volume of ferrous sulfate to titrate two blanks of potassium dichromate without soil (in mL)

**mLS** = the volume of ferrous sulfate to titrate the sample (in mL)

0.5 = the normality of the ferrous sulfate titrant

3. Substitution into the equation for OC results in the following.

$$\text{OC (\%)} = (\text{mLB} - \text{mLS}) \times 0.5 \times [0.003 \text{ (g.C.meq}^{-1})] \times \text{cf} \div \text{g soil} \times 100$$

Using a correction factor (cf) of 1.3 results in:

$$\text{OC (\%)} = (\text{mLB} - \text{mLS}) \times 0.195 \div \text{g soil}$$

4. Since there is approximately 58% C in soil organic matter, soil OM can be determined as:

$$\begin{aligned} \text{OM, \%} &= (\text{mLB} - \text{mLS}) \times 0.195 \div \text{g soil} \times (100 \text{ g OM} \div 58 \text{ g OC}) \\ &= (\text{mLB} - \text{mLS}) \times 0.336 \div \text{g soil} \end{aligned}$$

## Analytical Performance

### Range and Sensitivity

1. This method is appropriate for mineral soils with OM content less than 6%. Soils with more than 6% OM should be analysed using a quantity of soil lower than 1 g to avoid difficulties in determining the colour change during the titration and consumed all the potassium dichromate, which oxidises the OM.
2. Each volumetric reading is recorded to the nearest 0.1 mL. This titrant volume change corresponds to 0.03% OM when 1 g of soil is used.

### Precision and Accuracy

1. To ensure good precision, you must be able to view the reddish-brown colour change at the endpoint to stop the titration. The aid of a stirring bar and a well-lit stirring plate can improve the ability to view the end point.
2. Results from 151 North American Proficiency Testing samples from 1999 through 2008 show median interlaboratory precision of 10.6% ranging from 4.6 to 28.4% (median absolute deviation/median x 100) with median soil OM of 2.0% ranging from 0.5 to 8.1%.
3. Intralaboratory precision from the University of Florida Extension Soil Testing Laboratory on control samples, analysed 15 times, averaged  $1.09 \pm 0.074\%$  (standard deviation/average x 100) for soils with OM ranging from 0.85% to 1.14%.

### Interferences

1. Ferrous iron and chloride in soil can result in positive errors in OM as these constituents can result in a reduction of  $\text{K}_2\text{Cr}_2\text{O}_7$ . However, ferrous iron should not be a problem in soils thoroughly air-dried since this step ensures almost complete oxidation of ferrous iron to ferric iron. Chloride should not be a concern in well drained soils devoid of any recent fertiliser application. Options for treating soils with high chloride can be found in Nelson and Sommers (1996).



2. Manganese oxide ( $\text{MnO}_2$ ) can compete with  $\text{K}_2\text{Cr}_2\text{O}_7$  in oxidizing soil OM resulting in negative errors in OM. This is a rare occurrence and only occurs with freshly precipitated  $\text{MnO}_2$ .

### Interpretation

Interpretation of OM concentrations regarding soil fertility can be found in Magdoff *et al.* (1996).

### Effects of Storage

1. Air-dried soil can be stored indefinitely without affecting these measurements.
2. Potassium dichromate and ferrous sulfate solutions should be prepared fresh every day

### Safety and Disposal

Potassium dichromate is defined as a hazardous chemical by the Resource Conservation and Recovery Act due to toxicity (USEPA, 1980a). A laboratory generating more than 100 kg of material in a month is considered a hazardous waste generator that needs to follow hazardous waste disposal protocols defined by the US Environmental Protection Agency (USEPA, 1980b).

### Additional Reference

Magdoff, F.R., M.A. Tabatabai, and E.D. Hanlon. 1996. Soil Organic Matter: Analysis and Interpretation. Soil Sci. Spec. Pub. No. 46:21-31.

Nelson, D.W. and L.E. Sommers. 1996. Total carbon, organic carbon, and organic matter. p. 961-1010. In D.L Sparks *et al.* (eds.) Methods of soil analysis. Part 3. Chemical Methods. SSSA Book Series No. 5, SSSA and ASA, Madison, WI.

U.S. Environmental Protection Agency. 1980a. Title 40. Protection of Environment. Part 261—Identification and Listing of Hazardous Wastes. Revised as of 1 July 2004. U.S. Gov. Print. Office, Washington, DC.

U.S. Environmental Protection Agency. 1980b. Title 40. Protection of Environment. Part 262—Standards applicable to generators of hazardous wastes. Revised as of 1 July 2003. U.S. Gov. Print. Office, Washington, DC.

Walkley, A. and I.A. Black. 1934. An examination of the Degtjareff method for determining soil organic matter, and a proposed modification of the chromic acid titration method. Soil Sci. 37:29-38



### 3 Linked SOPs

SOP Code	SOP Function
EU_TRUE_SOP_025	Soil nutrient analyses

### 4 Disclaimer

The information presented here has been thoroughly researched and is believed to be accurate and correct. However, the authors cannot be held legally responsible for any errors. There are no warranties, expressed or implied, made with respect to the information provided. The authors will not be liable for any direct, indirect, special, incidental or consequential damages arising out of the use or inability to use the content of this publication.

### 5 Copyright

© All rights reserved. Reproduction and dissemination of material presented here for research, educational or other non-commercial purposes are authorised without any prior written permission from the copyright holders provided the source is fully acknowledged. Reproduction of material for sale or other commercial purposes is prohibited.

### 6 Citation

Please cite this report as follows:

Karavidas I., Ntatsi G., Savvas D. (2018). Standard Operating Procedure 26: Measuring Soil Organic Carbon (OC). Developed by the EU-H2020 project TRUE ('Transition paths to sustainable legume-based systems in Europe'), funded by the European Union's Horizon 2020 Research and Innovation programme under Grant Agreement Number 727973.