




TRansition paths to sUustainable
legume-based systems in EUrope

Isolation of Single Rhizobia Strains

Date: May, 31st 2018



This Project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 727973

www.true-project.eu



Content

1 Summary Information	3
1.1 Partner Summary	3
1.2 SOP Summary	4
2 Protocol Steps	5
3 Linked SOPs	7
4 Disclaimer	7
5 Copyright	7
6 Citation	7



1 Summary Information

1.1 Partner Summary

SOP Code	EU-TRUE_SOP_002
TRUE Partner Acronym	JHI
Primary Author	James, E.K. (Euan.James@hutton.ac.uk)
Other Authors	Maluk, M., Iannetta, P.P.M.
Linked Reference and Hyperlink (if available)	NA
Associated files to use with the SOP [and function]	NA



1.2 SOP Summary

Title

Isolation of Single Rhizobia Strains

Brief description

This details the isolation of single (pure) rhizobial strains from legume nodules using standard sterile techniques. Pure strains are essential for determining the symbionts involved in nodulation and Biological Nitrogen Fixation (BNF), especially if they are to be screened for their effectiveness at BNF and promoting plant growth with the ultimate aim of defining them as potential commercial inoculant strains.

2 Protocol Steps

Isolation of Single Rhizobia Strains

1. Prepare 70% ethanol, 2.5 % Na hypochlorite (NaClO; 2.5% active chlorine)¹, sterile distilled water. Sterile 1.5 mL tubes and sterile plastic pestles². Yeast mannitol agar (YMA³) plates.
2. For dry nodules, soak over-night in sterile distilled water.
For frozen nodules, allow them to thaw out.
3. Place nodules in 1.5 mL tube.
4. Surface sterilize nodules with 70 % ethanol for 1 min.
5. Sterilize nodules with 2.5 % Na hypochlorite for 3 min.
6. Wash with sterile water 3 times.
7. Squash the nodules in 1.5 mL tube using sterile plastic pestle.
8. Disperse the nodule juices across the YMA plate using sterile plastic pestle.
9. Use the same plastic pestle to spread onto next fresh YMA plate.
10. Allow the nodule juice to sink into the plate.
11. Incubate upside down (agar side upper-most), overnight at 28°C. Do not enclose/seal the plates in any way.⁴
12. Check for growth of colonies, and select individual white/pale and mucilaginous looking colonies for further growth. NB - avoid dark red colonies/growths, or those that look 'dry' or fungal. Aim to select one colony and streak this onto a fresh YMA plate. Again, incubate overnight at 28 °C.

¹ Household bleach is generally *ca.* 12.5 % chlorine, therefore dilute this x5.

² Pellet pestles Z359947-100EA Sigma-Aldrich;

<http://www.sigmaaldrich.com/catalog/product/sigma/z359947?lang=en®ion=GB>.

³ **YMA** - For 1 L: 10 g of mannitol; 0.5 g of glutamate; 0.5g of K₂HPO₄; 0.1 g of MgSO₄·7H₂O; 0.05 g of NaCl; 1 mL (40 g L⁻¹) CaCl₂; 1 mL (4 g L⁻¹) FeCl₃; 1g of yeast extract; 15 g technical grade agar. Add distilled water to 800 mL and pH to 6.8. Make volume to 1L. Autoclave, and whilst still warm (~60 °C), add 10 mL L⁻¹ of filter sterilised Congo Red Stock Solution (sterile; 0.25 % [w/v]). Dispense ~20 ml *per* 9 cm diameter Petri dish.

⁴ Fresh nodule extracts will yield cultures that generally grow well overnight. However, isolation from the other nodules types (dry or frozen), may take up to 3-4 days to appear as colonies that are large enough to discern their phenotype.



-
13. Once re-grown, take a single colony from the pure isolate on YMA plate and inoculate a sterile 5 ml **TY** broth⁵.
 14. Grow the TY at 28 °C (with shaking). After 24h (at log phase), take 0.9 mL of culture and combine with 0.9 mL of sterile 50% glycerol. Mix and freeze in liquid-N₂. Store the glycerol stock at -80 °C.
 15. The stock may now be used directly for bacterial DNA extractions and/or for PCR or infection tests.

⁵ **TY** - For 1 L: 5 g of tryptone; 3 g of yeast extract; 0.913 g of CaCl₂·2H₂O. Add distilled water to almost 1L, pH to 6.8, make volume to 1L. Autoclave. Dispense into 30 mL sterile universals ([Sterilin™ Polypropylene 30mL Universal Containers](#); Catalog # : 128A; ThermoFisher) as 5 mL aliquots.

3 Linked SOPs

SOP Code	SOP Function
EU_TRUE_SOP_001	Measuring ‘%Ndfa’ and ‘Biological Nitrogen Fixation’
EU_TRUE_SOP_003	Identification of root nodule isolates by 16S-rRNA PCR

4 Disclaimer

The information presented here has been thoroughly researched and is believed to be accurate and correct. However, the authors cannot be held legally responsible for any errors. There are no warranties, expressed or implied, made with respect to the information provided. The authors will not be liable for any direct, indirect, special, incidental or consequential damages arising out of the use or inability to use the content of this publication.

5 Copyright

© All rights reserved. Reproduction and dissemination of material presented here for research, educational or other non-commercial purposes are authorised without any prior written permission from the copyright holders provided the source is fully acknowledged. Reproduction of material for sale or other commercial purposes is prohibited.

6 Citation

Please cite this report as follows:

James, E.K., Maluk, M., Iannetta, P.P.M. (2018). *Isolation of Single Rhizobia Strains*. Developed by the EU-H2020 project TRUE (‘Transition paths to sustainable legume-based systems in Europe’), funded by the European Union’s Horizon 2020 Research and Innovation programme under Grant Agreement Number 727973. Available online at: www.true-project.eu.